

Come and see us next week at Dioxin 08!

Mass Spec UK will be exhibiting at the International Dioxin/POP's at the International Conference Centre in Birmingham from the 18th to 21st August.

On show will be our new Konik Pestilizer GC MS K2 system, which will be the first time the Konik has been on show in the UK.

So if you are attending Dioxin 08, we would love to see you. Please stop by Booth 1 for a chat and you can see the Konik Pestilizer.



Introducing another new Konik system.



In our last newsletter we informed you that we are now the UK distributor for Konik's range of analytical products. Konik's GC, HPLC, GC/MS and LC/MS systems are available with service and support directly from Mass Spec UK.


A new HPLC-GC-MS K2 Aqualizer system has just been introduced as another 'all-in-one' system in the growing range of Konik K2 systems. This system is primarily for the analysis of pesticides in water, but can also be configured for any other HPLS/MS or GC/MS application required.

Need an affordable benchtop GC-Tof instrument?

Many businesses are having to tighten the purse strings at the moment and if you need a quality benchtop GC-Tof instrument that is affordable, then your wait is over.

Mass Spec UK can now supply an SAI benchtop GC-Tof system that is ideal for those on a tight budget.

For more information on our affordable benchtop and to see if it will meet your requirements please contact us today on 0161 785 0828.



Service and Support now available for Thermo Finnigan LCQ

Mass Spec UK are pleased to announce that service and support for the Thermo Finnigan LCQ range of instruments is now available from our associate company, Labmend. Service contracts, PM's and paid call outs are available for LC Q Classic, Duo, Deca XP and Surveyor HPLC systems.

HINTS AND TIPS

Are you having trouble with your Mass Spec not pumping down after venting? If your rotary/roughing pump is fitted with a molecular sieve then the problem may be caused by excessive moisture/solvent vapour in the sieve.

To fix the problem simply bake out the pellets in the sieve or replace with new pellets.

To find out more contact Mass Spec UK on 0161 785 0828 today.

REFER A FRIEND/COLLEAGUE

Do you know someone who may find our newsletter interesting? You don't need to keep forwarding our newsletter on, they can receive one too.

Simply send an email to barrien@massspecuk.ltd.uk with 'Refer a friend' in the subject line and your friends details in the email, and we will make sure they receive the next issue of our newsletter.

sales@massspecuk.ltd.uk | www.massspecuk.ltd.uk | service@massspecuk.ltd.uk

Mass Spec UK Ltd, Highfield Street, Oldham, OL9 6DT | Tel: 0161 785 0828 | Fax: 0161 785 0838