

EXCITING NEW PARTNERSHIP WITH KONIK!

Mass Spec UK are very pleased to announce that we have been appointed as the official UK distributor for Konik's exciting range of scientific and analytical products.

From June 2008 Mass Spec UK will be responsible for the sales, service and support for the Konik range of GCs, LCs, HPLCs and Mass Spectrometer systems.

Konik is a leading European designer and manufacturer of complete and innovative solutions to meet customers' demanding requirements. The Konik Pestilizer is the first dedicated HPLCxHRGC-MS system (other dedicated systems to follow) and is available now!

The Konik Pestilizer is a UK first. Using the patented K2 interface, it allows for simplified sample preparation, minimum use of solvents, lower detection limits and a reduction in analysis costs. The use of the Konik multi-application Robokrom auto sampler allows for full automation of the Pestilizer HPLCxHRGC-MS system.

If you are not looking for a dedicated analyser, the modular and flexible range of the Konik's products allows us to offer instrumentation configured to your exact requirements. e.g: HPLC, LC, auto sampler can be added to a GC/MS system as and when required.



INTRODUCTORY SPECIAL OFFER!

OFFER 1

To celebrate our new partnership, Konik are offering a 'Fly, Try and Buy' programme.

Fly to the Konik application laboratories in Barcelona, and Konik will carry out a full day's demo at no cost!

OFFER 2

If you can't spare the time to fly to Barcelona, Mass Spec UK are offering a free feasibility study.

Just send us your samples, specify your detection limits and a thorough feasibility study will be carried out at no charge.

Contact us today on **0161 785 0828** for more information and how to apply for our special offers.

COME AND SEE US AND OUR NEW KONIK PESTILIZER

The 17th - 22nd of August are dates for your diary. We will be exhibiting at the Dioxin 2008 exhibition at the International Conference Centre in Birmingham and will be demonstrating the new Konik Pestilizer

This will be the first time the Konik pestilizer has been on show in the UK and is not an event to be missed.

More information on the event can be found at www.dioxin2008.org.

FAREWELL KEVIN



Some of you may be aware that we have recently said goodbye to Kevin Jones. Kevin set up Mass Spec UK with Barrie Nixon over 10 years ago.

We are sure that everyone who knows Kevin will join us in wishing him all the best for his future endeavors.

HINTS AND TIPS

Is your detector looking tired or at the end of its useful operating life? Looking for a better alternative to the original?

A replacement ETP can give an increase in sensitivity and performance for your Agilent 5973/5975 GC/MS system.

Replacement ETP detectors are available direct from Mass Spec UK. Please quote the following part numbers :

☞ MS689 ETP detector with mounting block. Mounting block is required for first-time-fit only

☞ MS690 ETP detector only

Mass Spec UK Ltd, Highfield Street, Oldham, OL9 6DT
Tel: 0161 785 0828 Fax: 0161 785 0838
sales@masspecuk.ltd.uk
service@massspecuk.ltd.uk
www.massspecuk.ltd.uk